



- Easy Neurac[®] testing and treatment
- Adaptable to all patients
- Full body suspension with three traverses

The Redcord Workstation Professional, with the triple traverse system and Redcord AXIS, enables the full range of Neurac[®] including full body weight-bearing positions.

The ONLY solution for Neurac[®] treatment

- ✓ Easy to get appropriate biomechanical positions with elastic cords
- ✓ Full range of suspension exercise progressions
- ✓ Quick adjustment of slings, straps and elastic cords
- ✓ Freedom of movement in the ropes
- ✓ Comfortable and effective management of pain
- ✓ Smooth gliding system for easy use



Perform simple self-mobilization exercises by releasing the ropes from the Redcord Trainer with the Redcord Rope Release.



Visit www.redcord.com to read about the Redcord Educational program.

Neurac treatment restores optimal and pain free movement patterns through neuromuscular stimulation.

Redcord® Workstation Professional

All you need for Neurac treatment



Ceiling Mounted Workstation

Art. nr.: 100003

For ceiling heights up to 350 cm

Dimensions:
Width: 78 cm
Length: 182 cm



Wall Mounted Workstation

Art. nr.: 100103

Dimensions:
Width: 201,6 cm
Depth: 228,5 cm
Height: 245 cm



Floor Standing Workstation

Art. nr.: 100203

Dimensions:
Width: 265,0 cm
Depth: 190,0 cm
Height: 2,52 cm - 2,56 cm (adjustable with feet)

Please contact your local distributor for advice about ceiling suspension options.

Workstation Professional content:

1 x Redcord Sliding Suspension System for ceiling, wall or floor mounting with 3 traverses.



3 x Redcord Trainer with straps



1 x Redcord AXIS with straps



2 x Redcord Wide Sling



2 x Redcord Narrow Sling



1 x Redcord Split Sling



1 pair Redcord PowerGrip



1 pair Redcord Elastic Cord 30 cm
1 pair Redcord Elastic Cord 60 cm
Black low resistance



1 pair Redcord Elastic Cord 30 cm
1 pair Redcord Elastic Cord 60 cm
Red high resistance



1 pair Redcord Rope, 30 cm
1 pair Redcord Rope, 60 cm



1 x Redcord Rope 5 m,
with hook



1 x Redcord Roll 15 x 50 cm
2 x Redcord Balance



Redcord Cord Lock
for Redcord AXIS



Redcord AXIS Ceiling
suspension with screws
3 x Carabiners



1 x Redcord Ceiling
Suspension, with screws



3 x Redcord Rope Clip
1 x Redcord Wall Caddy for slings, ropes



1 x Redcord Trainer/Mini
User Instructions
1 x Redcord Exercise Poster



1 x Redcord AXIS Exercise poster and
1 x Redcord AXIS User Manual



1 x Redcord Rope Release



Redcord Workstation is designed in compliance with the essential requirements for medical devices (Class 1, rule 1), and is labeled with the CE mark according to the manufacturer's declaration.